

OEC 9600[™]

Digital Mobile C-Arm



OEC 9600 Series Benefits

Video Imaging System

- 9' / 6' / 4.5' Imaging Intensifier - *New**
- 4.5' = 55 Lp / cm
- 6' = 48 Lp / cm
- 9' = 44 Lp / cm
- 32:1 Contrast Ratio
- 85% = QDE (typical)

X-Ray Tube - *New**

- Eureka - rotating anode
- 0.3mm - 1.0mm focal spots
- Focal spot kW rating: 0.3 - 5.3 kW, 1.0 - 41.0 kW
- Anode heat capacity: 300,000 HU
- Anode cooling rate: 70,000 HU per minute
- Housing heat capacity: 1,250,000 HU
- Anode diameter: 3"
- Anode angle: 12°

Collimator

- Two sets of opposing shutters
- Continuously variable
- Rotatable
- Remote control from the control panel

- **Great Images**
- **Low Dose**
- **Maneuverable**
- **Simple Operation**
- **Large Field of View**
- **Rapid Hard Copy**

